

Development of a reference method

to quantify the viral load

Yan Yan Beer^{1,2}, Andre Henrion¹, Melanie Brinkmann^{2,3}, Gavin O'Connor^{1,4}

¹Biochemistry, Physikalisch-Technische Bundesanstalt, 38116 Braunschweig, Germany

²Viral Immune Modulation Research Group, Helmholtz Centre for Infection Research, 38124 Braunschweig, Germany

³Institute of Genetics, Technische Universität Braunschweig, 38106 Braunschweig, Germany


⁴Institute of Biochemistry Biotechnology and Bioinformatics, Technische Universität Braunschweig, 38106 Braunschweig, Germany

BACKGROUND

Accurate and comparable measurements are essential to support the diagnosis of viral disease. Consistent measurement results can be achieved by making them SI-traceable through the development of reference measurements. DNA and Protein serve as measurands of interest.

ASSUMPTION


Viral Load / ml
 $\hat{=}$
 Viral Particle Amount / ml
 $\hat{=}$
 Viral Genome Copies / ml
 $\hat{=}$
 Viral Protein Quantities / ml



IMPACT

- SI-traceable method enables comparable measurements
- Alternative method to
 - diagnose viral disease and to determine viral load
 - identify false positive/negative results from already existing methods

PROJECT OUTLINE



Virus Preparation

- Virus production
- Virus purification
- Virus titration using qPCR

Proteomics Study

- Protein sample preparation
- LC-MS measurements
- Bioinformatically analysis to choose peptide of interest

Method Development

- Internal standard calibration
- Choose/test labeling method
- Assay verification
- Uncertainty measurements



References:

Arsene, Cristian G., et al. "Protein Quantification by Isotope Dilution Mass Spectrometry of Proteolytic Fragments : Cleavage Rate and Accuracy." *Analytical Chemistry*, vol. 80, no. 11, 2008, pp. 4154–60.
 Burkitt, William I., et al. "Toward Systeme International d'Unité-Traceable Protein Quantification: From Amino Acids to Proteins." *Analytical Biochemistry*, vol. 376, no. 2, 2008, pp. 242–51, doi:10.1016/j.ab.2008.02.010.
 Figures and icons were created with BioRender.com